

## In Situ Nick Translation (ISNT) for the Detection of Apoptotic Cells in Skin

1. Make cryo sections of human skin (7  $\mu$ m thick) and mount them on poly-L-lysine coated cover glasses. Let the sections dry for 30 min at room temperature.
2. Fixation in acetone (precooled) for 10 min at 4°C.
3. Wash in PBS for 5 min at r.t.
4. ISNT-reaction: 100  $\mu$ l/ cover glass; 40 min at r.t.: amount for 1 ml reaction mixture

3 $\mu$ M FITC-12-dUTP (Boehringer 1373242):	3.75 $\mu$ l (0.8 mM stock solution)
(or 3 $\mu$ M Biotin-16-dUTP)	
3 $\mu$ M dGTP	7.5 $\mu$ l (0.4 mM stock)
3 $\mu$ M dATP	7.5 $\mu$ l (0.4 mM stock)
3 $\mu$ M dCTP	7.5 $\mu$ l (0.4 mM stock)
DNA-Polymerase I (Boehringer 642711) (endonuclease-free)	50 units/ml (10 $\mu$ l of 5 u/ $\mu$ l stock)
10x reaction buffer (50 mM Tris/HCl pH7.5, 10 mM MgCl <sub>2</sub> , 0.1 mM DTT)	100 $\mu$ l
A. dest nuclease-free	ad 1 ml
5. Washing with PBS: 3 times for 5 min at r.t.
6. Protein block: 30 min at r.t. with 10% FCS in PBS
7. Incubation with peroxidase-conjugated anti-FITC (Boehringer; 1:25 in block solution; 30 min at 37°C)
8. Washing with PBS: 3 x 5 min at r.t.
9. Metal-enhanced DAB-staining (Pierce, 34065): 100  $\mu$ l/cover glass: Incubation about 5 - 20 min (r.t.)

Mounting and conventional microscopy (maybe after hemalaun counterstaining)

Alternative: (if Biotin-16-dUTP was used):
7. Incubation with FITC- or Texas Red conjugated streptavidin (Amersham)
8. Washing with PBS: 3 or 4 times for 5 min at r.t.
9. Mounting and fluorescence or laser scanning microscopy