syngo.share / Release VA25B / 2018-02-26 / Revision 2690

# syngo.share import

User Manual

Distributed by:

**SIEMENS**

# Language Directory

Dansk	3
Deutsch	14
Ελληνικά	24
English	34
Español	47
Français	57
Magyar	69
Italiano	81
Nederlands	91
Norsk (nynorsk)	101
Polski	112
Português	124
Türkçe	136

syngo.share / Release VA25B / 2018-02-26

# syngo.share import

Brugermanual

Distributed by:

**SIEMENS**

## Fordelene ved syngo.share

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ITH icoserve technology for healthcare GmbH  
Innrain 98  
A-6020 Innsbruck  
Austria

2018

**Distributed by**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
[siemens.com/healthcare](http://siemens.com/healthcare)

**Legal Manufacturer**

ITH icoserve technology for healthcare GmbH  
Innrain 98  
6020 Innsbruck  
Austria  
Phone: +43 512 89059

*syngo.share* / Release VA25B / 2018-02-26

# *syngo.share* import

User Manual

Distributed by:

**SIEMENS**

## Benefits of *syngo.share*

- Consistent, company-wide archive for all clinical pictures, multimedia, and documents.
- A flexible mechanism for the acquisition of all your clinical data (DICOM-/HL7-interfaces, import of data, video recordings, print interfaces, etc.).
- Uniform software for viewing the most common clinical data formats. Special framework conditions can be adjusted (display, lighting conditions) to facilitate diagnoses. According to the working environment, a bright or dark background can be enabled.
- Authorization system based on functions and positions
- Continue working while doing long-lasting tasks
- Integration with existing an HIS (Hospital Information System)

# Contents

- 1 About *syngo.share* import 37**
  - 1.1 Intended Use . . . . . 37
  - 1.2 System Prerequisites . . . . . 37
  
- 2 Safety Notices 38**
  - Data Import . . . . . 38
  - Medical Education . . . . . 38
  - Calculation Deviations . . . . . 38
  - Units . . . . . 38
  - Access Restriction . . . . . 39
  - Effects of External Programs . . . . . 39
  - Printing of Images . . . . . 39
  
- 3 First Steps 39**
  - 3.1 User Security . . . . . 39
  - 3.2 Keystroke Combination . . . . . 40
  - 3.3 Search for Patients and Documents . . . . . 40
    - 3.3.1 Advanced Search with Wildcards . . . . . 40
  - 3.4 Document Search and Results (*syngo.share* view, *syngo.share* webview) . . . . . 40
  - 3.5 Keyword Tagging . . . . . 41
  - 3.6 Export/local Saving . . . . . 41
    - 3.6.1 Exporting with *syngo.share* view . . . . . 41
    - 3.6.2 Exporting with *syngo.share* webview . . . . . 41
  - 3.7 Control State of Actions . . . . . 42
  
- 4 Application *syngo.share* import 42**
  - 4.1 Functionality of *syngo.share* import . . . . . 42
  
- EC Declaration of Conformity *syngo.share* import 43**
  
- Australian Conformity *syngo.share* import 44**

# 1 About *syngo.share* import

An application to import and tag data from various sources such as file systems, CD/DVD-ROM, USB drives, DICOMDIR CDs, digital video cameras, digital photography cameras, document and DICOM scanners, and DICOM stores to *syngo.share* VNA or an alternate DICOM archive.

## 1.1 Intended Use

*syngo.share* import is software for the electronic display, conversion and transfer of digital data, without controlling or altering the functions of any connected medical devices.

## 1.2 System Prerequisites

The following operating systems are supported:

- Windows 10 Professional / Enterprise Edition 64 Bit
- Windows 10 Professional / Enterprise Edition 32 Bit
- Windows 8.1 Professional / Enterprise Edition 64 Bit
- Windows 8.1 Professional / Enterprise Edition 32 Bit
- Windows 7 Professional / Enterprise Edition 64 Bit SP1

The minimum resolution for *syngo.share* import is 1280x1024.

## 2 Safety Notices

### Data Import

*syngo.share* import is only intended for the importing of data from various sources.

**Do not use *syngo.share* import for direct diagnosis.**

### Medical Education

The use of *syngo.share* import should be granted to qualified staff members only.

### Calculation Deviations

Any calculations and calibrations not made in *syngo.share* import will have a level of precision that correlates to the precision of the program of that in which they were created—not of *syngo.share* import. Additionally, the level of precision is only as precise as the data collected by the modality used. This means that each time a modality is connected, it should be calibrated. The calibration of the device and interpretation of calculations used in a diagnosis must be reviewed by qualified staff members for plausibility. Deviations may be caused by the image itself—in particular for conventional X-ray scans (CR) e.g. pelvis or hips. The representation of the distance between the part of the body being recorded and the imaging plate could deform the length on a two-dimensional X-ray image representation. Even if there is a point of reference included in the scan (e.g. a ball), there will always be a deviation if the point of reference is not positioned precisely at the height of the body part being calculated at the time of the scan. The area of polygons whose edges overlap will not be calculated correctly in each case. In addition, the extrapolation algorithm used in *syngo.share* import for displaying images can occasionally lead to imprecise rendering. Calculations of such images could be inaccurate.

The following defines how calculations will be rounded:

- Length and surface area: to the nearest tenth
- Density (HU): to the nearest whole number

### Units

The application *syngo.share* import conforms to Directive 80/181/EEC, which contains specifications related to units of calculation. Calibrated lengths are to be calculated using millimeters. Surface areas are to

be calculated in square millimeters (mm<sup>2</sup>). Plane angles are to be calculated in degrees. Densities will be calculated using the industry-standard unit, Hounsfield units (HU), which are not regulated by Directive 80/181/EEC.

### Access Restriction

Be sure to log out of *syngo.share* import when finished working in order to avoid unauthorized access to the system.

### Effects of External Programs

The software manufacturer does not ensure that third-party applications will not negatively affect the system. Any external programs accessed by *syngo.share* import have standards which correspond to their respective producers.

### Printing of Images

When printing images using the image caption option “As camera”, it is possible that some labels will be hidden from view due to the physically limited space on the page. As a result, possible transformations such as mirroring, rotation, and applied presentation states etc.; and even patient information may not be visible on the printout. If a series has active presentation states in both the view and print layouts, deactivation of the presentation states in the view layout will also deactivate it in the print layout of the series.

## 3 First Steps

This chapter contains the most important features of *syngo.share* as well as operational instructions. The applications *syngo.share* view, *syngo.share* print, and *syngo.share* import can be started from the Windows® start menu. The web applications *syngo.share* webview and *webbox* run inside a web browser at (e.g. <https://syngo.share/>).

### 3.1 User Security

A valid username and password is required by *syngo.share*, and will be assigned by a system administrator at your institution. The allocated roles and permissions determine which features will be enabled for a

given user. An inactive application will be automatically locked to ensure that there is no unauthorized access. To continue working where you left off, simply log in again and the program will resume at the point just before the automatic lock was activated.

## 3.2 Keystroke Combination

You can use keystroke combinations to initiate functions. The shortcuts are shown in the context menu or in the tool tip of the particular function. A complete list of all keystroke combinations can be found in the online help manual.

## 3.3 Search for Patients and Documents

To begin a search, insert the appropriate search criteria into the search field and press the "Search" button. If your search criteria is insufficient, the "Search" button will remain inactive. To remove all search criteria quickly, press the "Reset" button. The fields of the search mask can be defined at the applications view and *syngo.share* webview by using the button "Configure search fields" and by (de)activating the needed search fields. This setting will be saved by the programs. The use of capital and small initial letters won't be regarded.

### 3.3.1 Advanced Search with Wildcards

The following characters are used to perform advanced searches with wildcards in the search interface:

- The sign "%" is a variable for one or more characters
- The underscore "\_" is a variable for exactly one optional character

#### Examples:

- "m%er" would result in: "Meier", "Mayer", "Moser", "Mariacher", etc.
- "m\_ier" would result in: "Meier", "Maier", etc.

## 3.4 Document Search and Results (*syngo.share* view, *syngo.share* webview)

The search results columns can be arbitrarily ordered by simply dragging them with the mouse pointer. The order you chose will be saved in your user profile and used in future searches.

To sort a column, click on the column header.

- In the web application, you can click on the arrow in the column header to show or hide columns. These settings can be configured in the context menu of the application *syngo.share* view. Your settings will be saved and available the next time you launch the program.

- The search results can be exported as a CSV file in the web applications *syngo.share* webview and *webadmin*.

To open a search result, double click it in the list of results:

- In the application *syngo.share* view, the patient record will be loaded into the document tree.
- In the applications *syngo.share* webview and *webbox*, the selected item will be shown in the viewer.

## 3.5 Keyword Tagging

In the applications *syngo.share* view, *syngo.share* print, *syngo.share* import, and *syngo.share* webview you can annotate documents. Special authorization from an administrator is needed for this functionality. Tags add another layer of description to documents such that they can be easily found in the future. A configuration determines whether the tags can be chosen from a catalog, dynamically added to a catalog, or entered as free text (e.g. comments). The following DICOM levels support keyword tags: studies, series, and images. Generic containers also support tagging. While searching using keywords, the results are shown independently of the level in which the keyword tag was added.

For further details for tagging read the online help.

## 3.6 Export/local Saving

In the applications *syngo.share* view, *syngo.share* webview, and *webbox*, there are export functions that allow the data to be saved locally. To export data, two different levels of permission are possible:

- Anonymous-only export
- Full permissions export

### 3.6.1 Exporting with *syngo.share* view

Exporting of studies and series in the document tree:

- Click "Export all series" button to export every file.
- Select the desired items and export them by using the button "Export selected series".

Exporting of individual images in the print area:

1. Drag the requested element out of the document tree and into the print area.
2. Click the "Export..." button.

### 3.6.2 Exporting with *syngo.share* webview

Exporting a selection in *syngo.share* webview:

1. Click on either "Export..." or "Export as..." in the context menu.
2. The export wizard will open where you can configure your export.

## 3.7 Control State of Actions

All desktop applications have a status area, which contain

- If the requested action was successful the state will show a green check mark.
- An exclamation point indicates that the action could not be completed successfully. Click on the status for detailed information.

# 4 Application *syngo.share* import

In the application *syngo.share* import, you can archive any kind of file (DICOM objects, pictures, videos, documents, etc). To use it, you must do the following steps:

- Create a patient context.
- Select the required data and import it.
- Select a save location.
- Click on the "Archive" button.

The status area will indicate whether the archival process was successful.

## 4.1 Functionality of *syngo.share* import

- In the tab "Filesystem", you can import data from disk drives and USB flash drives:
  - Double click the desired files.
  - Drag them into the document menu bar.
- You can drag files directly from Windows® Explorer.
- To import patient CDs, choose the "Patient CD" tab.
- If there is a scanner connected, you can digitalize analogue data such as clinical reports or printed X-ray images. They can be edited after scanning and will be added to the archive.
- If there is a video source connected, you can record video streams in the "Video" tab. Additionally, it is possible to take a screenshot during a playing videostream by pressing the **F5** key.
- You can view a document in the document list with the viewer by double clicking it or by using the actions menu.
- You assign keyword tags for documents in the document list.

## EC Declaration of Conformity

### Product Details

Product Name	syngo.share import
Type	VA25B
Model Number	11250002
Classification according to Directive 93/42/EEC, Annex IX	Class I

### Manufacturer

Address	ITH icoserve technology for healthcare GmbH Innrain 98 6020 Innsbruck Austria Phone: +43 512 89059
---------	--

### Assessment Details

Procedure according to Directive 93/42/EEC	Annex VII
--	-----------

### Applied Standards

EN ISO 13485:2012 / AC:2012  
EN ISO 9001:2008/AC:2009  
EN ISO 14971:2012  
EN 62304:2006/AC:2008  
EN 62366:2008

*We declare under sole responsibility that the products described above are in compliance with Directive 93/42/EEC, if their installation and use occurs corresponding to the user and product documentation. This declaration supersedes any declaration previously issued for this product.*

Innsbruck, 2018-02-26

Dr. Bernhard Hirsch  
Managing Director

Mag. (FH) Herbert Plankl  
Head of Quality & Documentation

EC Declaration of Conformity

**Manufacturer's Declaration of Conformity**  
Australian Therapeutic Goods (Medical Devices) Regulations 2002  
Full Quality Assurance Procedure

This is a declaration made in accordance with the requirements of Clause 1.8 of Schedule 3 of the **Australian Therapeutic Goods (Medical Devices) Regulations 2002** relating to the stated device.

**Reference**

syngo.share import

**Legal Manufacturer Business Address**

ITH icoserve technology for healthcare GmbH  
Innrain 98  
6020 Innsbruck  
Austria  
Phone: +43 512 89059

**Facility**

ITH icoserve technology for healthcare GmbH  
Innrain 98  
6020 Innsbruck  
Austria  
Phone: +43 512 89059

**Medical Device(s)**

syngo.share import VA25B  
11250002

**Authorised Signatory**

2018-02-26

2018-02-26

---

Dr. Bernhard Hirsch  
Managing Director

---

Mag. (FH) Herbert Plankl  
Head of Quality & Documentation

**Declaration of Conformity AU**

- ☐ Classification (per TGA)
  - Class I
  - (Schedule 2, Part 4, Rule 4.1 of the Therapeutic Goods (Medical Devices) Regulations 2002)
- ☐ GMDN Code and Term
  - 41670: Radiology picture archiving and communication system application software
- ☐ Scope of Application
  - This declaration of conformity covers *syngo.share* import:
  - Model Number: 11250002
  - Product Description: *syngo.share* import VA25B

Each kind of medical device to which the Full Quality Assurance Procedures have been applied complies with the applicable provisions of the essential principles, the classification rules, at each stage, from the design of the device until its final inspection before being supplied.

### **Full Quality Management System Certificate**

EN ISO 13485:2012  
Certificate No.: 17062224

### **Standards Applied**

(European Harmonised Standard)

- EN ISO 9001:2008, AC:2009
- EN ISO 13485:2012, AC:2012
- EN ISO 14971:2012
- EN 62304:2006, AC:2008
- EN 62366:2008

ITH icoserve technology for healthcare GmbH  
Innrain 98  
A-6020 Innsbruck  
Austria

2018

**Distributed by**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
[siemens.com/healthcare](http://siemens.com/healthcare)

**Legal Manufacturer**

ITH icoserve technology for healthcare GmbH  
Innrain 98  
6020 Innsbruck  
Austria  
Phone: +43 512 89059

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