syngo.share / Release VA25B / 2018-02-26 / Revision 2690

syngo.share print

User Manual

Distributed by:

SIEMENS

Language Directory

Dansk	3
Deutsch	14
Ελληνικά	23
English	33
Español	44
Français	53
Magyar	65
Italiano	76
Nederlands	85
Norsk (nynorsk)	94
Polski	105
Português	116
Türkçe	127

syngo.share / Release VA25B / 2018-02-26

syngo.share print

Brugermanual

Distributed by:

SIEMENS

Fordelene ved *syngo.share*

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ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2018

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

syngo.share / Release VA25B / 2018-02-26

syngo.share print

User Manual

Distributed by:

SIEMENS

Benefits of *syngo.share*

- Consistent, company-wide archive for all clinical pictures, multimedia, and documents.
- A flexible mechanism for the acquisition of all your clinical data (DICOM-/HL7-interfaces, import of data, video recordings, print interfaces, etc.).
- Uniform software for viewing the most common clinical data formats. Special framework conditions can be adjusted (display, lighting conditions) to facilitate diagnoses. According to the working environment, a bright or dark background can be enabled.
- Authorization system based on functions and positions
- Continue working while doing long-lasting tasks
- Integration with existing an HIS (Hospital Information System)

Contents

- 1 About *syngo.share* print** **34**
- 1.1 Intended Use 34
- 1.2 System Prerequisites 34

- 2 Safety Notices** **35**
- Data Import 35
- Medical Education 35
- Calculation Deviations 35
- Units 35
- Access Restriction 36
- Effects of External Programs 36
- Printing of Images 36

- 3 First Steps** **36**
- 3.1 User Security 36
- 3.2 Keystroke Combination 37
- 3.3 Search for Patients and Documents 37
- 3.3.1 Advanced Search with Wildcards 37
- 3.4 Document Search and Results (*syngo.share* view, *syngo.share* webview) 37
- 3.5 Keyword Tagging 38
- 3.6 Export/local Saving 38
- 3.6.1 Exporting with *syngo.share* view 38
- 3.6.2 Exporting with *syngo.share* webview 38
- 3.7 Control State of Actions 39

- 4 Application *syngo.share* print** **39**

- EC Declaration of Conformity *syngo.share* print** **40**

- Australian Conformity *syngo.share* print** **41**

1 About *syngo.share print*

syngo.share print is an application that provides a virtual printer such that an arbitrary application can print directly to *syngo.share VNA* or an alternate DICOM archive as a PDF/A or a DICOM encapsulated PDF.

1.1 Intended Use

syngo.share print is software for the electronic conversion and transfer of digital data, without controlling or altering the functions of any connected medical devices.

1.2 System Prerequisites

syngo.share print can be used with the following operating systems:

- Windows 10 Professional / Enterprise Edition 64 Bit
- Windows 10 Professional / Enterprise Edition 32 Bit
- Windows 8.1 Professional / Enterprise Edition 64 Bit
- Windows 8.1 Professional / Enterprise Edition 32 Bit
- Windows 7 Professional / Enterprise Edition 64 Bit SP1

syngo.share print needs a resolution of at least 1280x1024.

2 Safety Notices

Data Import

syngo.share print is only intended for the importing of data from various sources.

Do not use *syngo.share print* for direct diagnosis.

Medical Education

The use of *syngo.share print* should be granted to qualified staff members only.

Calculation Deviations

Any calculations and calibrations not made in *syngo.share print* will have a level of precision that correlates to the precision of the program of that in which they were created—not of *syngo.share print*. Additionally, the level of precision is only as precise as the data collected by the modality used. This means that each time a modality is connected, it should be calibrated. The calibration of the device and interpretation of calculations used in a diagnosis must be reviewed by qualified staff members for plausibility. Deviations may be caused by the image itself—in particular for conventional X-ray scans (CR) e.g. pelvis or hips. The representation of the distance between the part of the body being recorded and the imaging plate could deform the length on a two-dimensional X-ray image representation. Even if there is a point of reference included in the scan (e.g. a ball), there will always be a deviation if the point of reference is not positioned precisely at the height of the body part being calculated at the time of the scan. The area of polygons whose edges overlap will not be calculated correctly in each case. In addition, the extrapolation algorithm used in *syngo.share print* for displaying images can occasionally lead to imprecise rendering. Calculations of such images could be inaccurate.

The following defines how calculations will be rounded:

- Length and surface area: to the nearest tenth
- Density (HU): to the nearest whole number

Units

The application *syngo.share print* conforms to Directive 80/181/EEC, which contains specifications related to units of calculation. Calibrated lengths are to be calculated using millimeters. Surface areas are to be

calculated in square millimeters (mm²). Plane angles are to be calculated in degrees. Densities will be calculated using the industry-standard unit, Hounsfield units (HU), which are not regulated by Directive 80/181/EEC.

Access Restriction

Be sure to log out of *syngo.share print* when finished working in order to avoid unauthorized access to the system.

Effects of External Programs

The software manufacturer does not ensure that third-party applications will not negatively affect the system. Any external programs accessed by *syngo.share print* have standards which correspond to their respective producers.

Printing of Images

When printing images using the image caption option “As camera”, it is possible that some labels will be hidden from view due to the physically limited space on the page. As a result, possible transformations such as mirroring, rotation, and applied presentation states etc.; and even patient information may not be visible on the printout. If a series has active presentation states in both the view and print layouts, deactivation of the presentation states in the view layout will also deactivate it in the print layout of the series.

3 First Steps

This chapter contains the most important features of *syngo.share* as well as operational instructions. The applications *syngo.share view*, *syngo.share print*, and *syngo.share import* can be started from the Windows® start menu. The web applications *syngo.share webview* and *webbox* run inside a web browser at (e.g. <https://syngo.share/>).

3.1 User Security

A valid username and password is required by *syngo.share*, and will be assigned by a system administrator at your institution. The allocated roles and permissions determine which features will be enabled for a

given user. An inactive application will be automatically locked to ensure that there is no unauthorized access. To continue working where you left off, simply log in again and the program will resume at the point just before the automatic lock was activated.

3.2 Keystroke Combination

You can use keystroke combinations to initiate functions. The shortcuts are shown in the context menu or in the tool tip of the particular function. A complete list of all keystroke combinations can be found in the online help manual.

3.3 Search for Patients and Documents

To begin a search, insert the appropriate search criteria into the search field and press the “Search” button. If your search criteria is insufficient, the “Search” button will remain inactive. To remove all search criteria quickly, press the “Reset” button. The fields of the search mask can be defined at the applications *syngo.share* view and *syngo.share* webview by using the button “Configure search fields” and by (de)activating the needed search fields. This setting will be saved by the programs. The use of capital and small initial letters won’t be regarded.

3.3.1 Advanced Search with Wildcards

The following characters are used to perform advanced searches with wildcards in the search interface:

- The sign “%” is a variable for one or more characters
- The underscore “_” is a variable for exactly one optional character

Examples:

- “m%er” would result in: “Meier”, “Mayer”, “Moser”, “Mariacher”, etc.
- “m_ier” would result in: “Meier”, “Maier”, etc.

3.4 Document Search and Results (*syngo.share* view, *syngo.share* webview)

The search results columns can be arbitrarily ordered by simply dragging them with the mouse pointer. The order you chose will be saved in your user profile and used in future searches.

To sort a column, click on the column header.

- In the web application, you can click on the arrow in the column header to show or hide columns. These settings can be configured in the context menu of the application *syngo.share* view. Your settings will be saved and available the next time you launch the program.

- The search results can be exported as a CSV file in the web applications *syngo.share* webview and *webadmin*.

To open a search result, double click it in the list of results:

- In the application *syngo.share* view, the patient record will be loaded into the document tree.
- In the applications *syngo.share* webview and *webbox*, the selected item will be shown in the viewer.

3.5 Keyword Tagging

In the applications *syngo.share* view, *syngo.share* print, *syngo.share* import, and *syngo.share* webview you can annotate documents. Special authorization from an administrator is needed for this functionality. Tags add another layer of description to documents such that they can be easily found in the future. A configuration determines whether the tags can be chosen from a catalog, dynamically added to a catalog, or entered as free text (e.g. comments). The following DICOM levels support keyword tags: studies, series, and images. Generic containers also support tagging. While searching using keywords, the results are shown independently of the level in which the keyword tag was added.

For further details for tagging read the online help.

3.6 Export/local Saving

In the applications *syngo.share* view, *syngo.share* webview, and *webbox*, there are export functions that allow the data to be saved locally. To export data, two different levels of permission are possible:

- Anonymous-only export
- Full permissions export

3.6.1 Exporting with *syngo.share* view

Exporting of studies and series in the document tree:

- Click "Export all series" button to export every file.
- Select the desired items and export them by using the button "Export selected series".

Exporting of individual images in the print area:

1. Drag the requested element out of the document tree and into the print area.
2. Click the "Export..." button.

3.6.2 Exporting with *syngo.share* webview

Exporting a selection in *syngo.share* webview:

1. Click on either "Export..." or "Export as..." in the context menu.
2. The export wizard will open where you can configure your export.

3.7 Control State of Actions

All desktop applications have a status area, which contain

- If the requested action was successful the state will show a green check mark.
- An exclamation point indicates that the action could not be completed successfully. Click on the status for detailed information.

4 Application *syngo.share print*

The application *syngo.share print* is a virtual printer that can quickly export your files directly into the archive. Any program that has a print function can use it. The original file is converted into a PDF file or embedded within a DICOM file. You can use *syngo.share print* with the following steps:

- Invoke the print function in whatever application you are using.
- Choose the printer "icoprint".
- Define a patient context in the new window:
 - Select a patient context from the worklist
 - Search for a patient
 - Providing patient data
- Define a save location.
- Click the "Archive" button to save the document permanently.

EC Declaration of Conformity

Product Details

Product Name	syngo.share print
Type	VA25B
Model Number	11250003
Classification according to Directive 93/42/EEC, Annex IX	Class I

Manufacturer

Address	ITH icoserve technology for healthcare GmbH Innrain 98 6020 Innsbruck Austria Phone: +43 512 89059
---------	--

Assessment Details

Procedure according to Directive 93/42/EEC	Annex VII
--	-----------

Applied Standards

EN ISO 13485:2012 / AC:2012
EN ISO 9001:2008/AC:2009
EN ISO 14971:2012
EN 62304:2006/AC:2008
EN 62366:2008

We declare under sole responsibility that the products described above are in compliance with Directive 93/42/EEC, if their installation and use occurs corresponding to the user and product documentation. This declaration supersedes any declaration previously issued for this product.

Innsbruck, 2018-02-26

Dr. Bernhard Hirsch
Managing Director

Mag. (FH) Herbert Plankl
Head of Quality & Documentation

EC Declaration of Conformity

Manufacturer's Declaration of Conformity
Australian Therapeutic Goods (Medical Devices) Regulations 2002
Full Quality Assurance Procedure

This is a declaration made in accordance with the requirements of Clause 1.8 of Schedule 3 of the **Australian Therapeutic Goods (Medical Devices) Regulations 2002** relating to the stated device.

Reference

syngo.share print

Legal Manufacturer Business Address

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

Facility

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

Medical Device(s)

syngo.share print VA25B
11250003

Authorised Signatory

2018-02-26

2018-02-26

Dr. Bernhard Hirsch
Managing Director

Mag. (FH) Herbert Plankl
Head of Quality & Documentation

Declaration of Conformity AU

- ☐ Classification (per TGA)
 - Class I
 - (Schedule 2, Part 4, Rule 4.1 of the Therapeutic Goods (Medical Devices) Regulations 2002)
- ☐ GMDN Code and Term
 - 41670: Radiology picture archiving and communication system application software
- ☐ Scope of Application
 - This declaration of conformity covers *syngo.share* print:
 - Model Number: 11250003
 - Product Description: *syngo.share* print VA25B

Each kind of medical device to which the Full Quality Assurance Procedures have been applied complies with the applicable provisions of the essential principles, the classification rules, at each stage, from the design of the device until its final inspection before being supplied.

Full Quality Management System Certificate

EN ISO 13485:2012
Certificate No.: 17062224

Standards Applied

(European Harmonised Standard)

- EN ISO 9001:2008, AC:2009
- EN ISO 13485:2012, AC:2012
- EN ISO 14971:2012
- EN 62304:2006, AC:2008
- EN 62366:2008

ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2018

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

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