syngo.share / Release VA25B / 2018-02-26 / Revision 2690

syngo.share view

User Manual

Distributed by:

SIEMENS

Language Directory

| | |
|-----------------|-----|
| Dansk | 3 |
| Deutsch | 17 |
| Ελληνικά | 30 |
| English | 43 |
| Español | 59 |
| Français | 71 |
| Magyar | 83 |
| Italiano | 95 |
| Nederlands | 107 |
| Norsk (nynorsk) | 119 |
| Polski | 133 |
| Português | 145 |
| Türkçe | 157 |

syngo.share / Release VA25B / 2018-02-26

syngo.share view

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Distributed by:

SIEMENS

Fordelene ved syngo.share view

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ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2018

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

syngo.share / Release VA25B / 2018-02-26

syngo.share view

User Manual English

Distributed by:

SIEMENS

Benefits of *syngo.share* view

- Consistent, company-wide archive for all clinical pictures, multimedia, and documents.
- A flexible mechanism for the acquisition of all your clinical data (DICOM-/HL7-interfaces, import of data, video recordings, print interfaces, etc.).
- Uniform software for viewing the most common clinical data formats. Special framework conditions can be adjusted (display, lighting conditions) to facilitate diagnoses. According to the working environment, a bright or dark background can be enabled.
- Authorization system based on functions and positions
- Continue working while doing long-lasting tasks
- Integration with existing an HIS (Hospital Information System)

Contents

| | | |
|----------|--|-----------|
| 1 | About <i>syngo.share</i> view | 46 |
| 1.1 | Intended Use | 46 |
| 1.2 | System Prerequisites | 46 |
| 2 | Safety Notices | 47 |
| | Image and Result Distribution | 47 |
| | Medical Education | 47 |
| | Calculation Deviations | 47 |
| | Units | 47 |
| | Access Restriction | 48 |
| | Effects of External Programs | 48 |
| | Printing of Images | 48 |
| 3 | First Steps | 49 |
| 3.1 | User Security | 49 |
| 3.2 | Keystroke Combination | 49 |
| 3.3 | Search for Patients and Documents | 49 |
| | 3.3.1 Advanced Search with Wildcards | 49 |
| 3.4 | Document Search and Results (<i>syngo.share</i> view, <i>syngo.share</i> webview) | 50 |
| 3.5 | Keyword Tagging | 50 |
| 3.6 | Export/local Saving | 50 |
| | 3.6.1 Exporting with <i>syngo.share</i> view | 51 |
| | 3.6.2 Exporting with <i>syngo.share</i> webview | 51 |
| 3.7 | Control State of Actions | 51 |
| 4 | Application <i>syngo.share</i> view | 51 |
| 4.1 | Display Mode | 51 |
| 4.2 | Load Area | 52 |
| 4.3 | View Area | 53 |
| | 4.3.1 Single-frame view | 53 |
| | 4.3.2 Movie mode | 53 |
| | 4.3.3 Metadata | 53 |
| | 4.3.4 Create and Manage collection folders | 54 |
| 4.4 | Print Area | 54 |
| | EC Declaration of Conformity <i>syngo.share</i> view | 55 |
| | Australian Conformity <i>syngo.share</i> view | 56 |

1 About *syngo.share* view

syngo.share view is a feature rich, lightweight desktop application to be used in clinical image and result distribution. It provides access to content from the *syngo.share* VNA, DICOM archives, XDS/I Affinity Domains, or a local file system. It integrates with EMR and clinical information systems to provide basic review and annotation support.

1.1 Intended Use

syngo.share view is software for the electronic display of data, without controlling or altering the functions of any connected medical devices.

1.2 System Prerequisites

The following operating systems are supported:

- Windows 10 Professional / Enterprise Edition 64 Bit
- Windows 10 Professional / Enterprise Edition 32 Bit
- Windows 8.1 Professional / Enterprise Edition 64 Bit
- Windows 8.1 Professional / Enterprise Edition 32 Bit
- Windows 7 Professional / Enterprise Edition 64 Bit SP1

The minimum resolution for *syngo.share* view is 1280x1024.

2 Safety Notices

Image and Result Distribution

syngo.share view is only intended for clinical image and result distribution.

Do not use *syngo.share* view for direct diagnoses.

Medical Education

The use of *syngo.share* view should be granted to qualified staff members only.

Calculation Deviations

Any calculations and calibrations not made in *syngo.share* view will have a level of precision that correlates to the precision of the program of that in which they were created—not of *syngo.share* view. Additionally, the level of precision is only as precise as the data collected by the modality used. This means that each time a modality is connected, it should be calibrated. The calibration of the device and interpretation of calculations used in a diagnosis must be reviewed by qualified staff members for plausibility. Deviations may be caused by the image itself—in particular for conventional X-ray scans (CR) e.g. pelvis or hips. The representation of the distance between the part of the body being recorded and the imaging plate could deform the length on a two-dimensional X-ray image representation. Even if there is a point of reference included in the scan (e.g. a ball), there will always be a deviation if the point of reference is not positioned precisely at the height of the body part being calculated at the time of the scan. The area of polygons whose edges overlap will not be calculated correctly in each case. In addition, the extrapolation algorithm used in *syngo.share* view for displaying images can occasionally lead to imprecise rendering. Calculations of such images could be inaccurate.

The following defines how calculations will be rounded:

- Length and surface area: to the nearest tenth
- Density (HU): to the nearest whole number

Units

The application *syngo.share* view conforms to Directive 80/181/EEC, which contains specifications related to units of calculation. Calibrated lengths are to be calculated using millimeters. Surface areas are to be

calculated in square millimeters (mm²). Plane angles are to be calculated in degrees. Densities will be calculated using the industry-standard unit, Hounsfield units (HU), which are not regulated by Directive 80/181/EEC.

Access Restriction

Be sure to log out of *syngo.share* view when finished working in order to avoid unauthorized access to the system.

Effects of External Programs

The software manufacturer does not ensure that third-party applications will not negatively affect the system. Any external programs accessed by *syngo.share* view have standards which correspond to their respective producers.

Printing of Images

When printing images using the image caption option “As camera”, it is possible that some labels will be hidden from view due to the physically limited space on the page. As a result, possible transformations such as mirroring, rotation, and applied presentation states etc.; and even patient information may not be visible on the printout. If a series has active presentation states in both the view and print layouts, deactivation of the presentation states in the view layout will also deactivate it in the print layout of the series.

3 First Steps

This chapter contains the most important features of *syngo.share* view as well as operational instructions. The applications *syngo.share view*, *syngo.share print*, and *syngo.share import* can be started from the Windows® start menu. The web applications *syngo.share webview* and *webbox* run inside a web browser at (e.g. <https://syngo.share/>).

3.1 User Security

A valid username and password is required by *syngo.share view*, and will be assigned by a system administrator at your institution. The allocated roles and permissions determine which features will be enabled for a given user. An inactive application will be automatically locked to ensure that there is no unauthorized access. To continue working where you left off, simply log in again and the program will resume at the point just before the automatic lock was activated.

3.2 Keystroke Combination

You can use keystroke combinations to initiate functions. The shortcuts are shown in the context menu or in the tool tip of the particular function. A complete list of all keystroke combinations can be found in the online help manual.

3.3 Search for Patients and Documents

To begin a search, insert the appropriate search criteria into the search field and press the “Search” button. If your search criteria is insufficient, the “Search” button will remain inactive. To remove all search criteria quickly, press the “Reset” button. The fields of the search mask can be defined at the applications view and *syngo.share webview* by using the button “Configure search fields” and by (de)activating the needed search fields. This setting will be saved by the programs. The use of capital and small initial letters won’t be regarded.

3.3.1 Advanced Search with Wildcards

The following characters are used to perform advanced searches with wildcards in the search interface:

- The sign “%” is a variable for one or more characters
- The underscore “_” is a variable for exactly one optional character

Examples:

- "m%er" would result in: "Meier", "Mayer", "Moser", "Mariacher", etc.
- "m_ier" would result in: "Meier", "Maier", etc.

3.4 Document Search and Results (*syngo.share view*, *syngo.share webview*)

The search results columns can be arbitrarily ordered by simply dragging them with the mouse pointer. The order you chose will be saved in your user profile and used in future searches.

To sort a column, click on the column header.

- In the web application, you can click on the arrow in the column header to show or hide columns. These settings can be configured in the context menu of the application *syngo.share view*. Your settings will be saved and available the next time you launch the program.
- The search results can be exported as a CSV file in the web applications *syngo.share webview* and *webadmin*.

To open a search result, double click it in the list of results:

- In the application *syngo.share view*, the patient record will be loaded into the document tree.
- In the applications *syngo.share webview* and *webbox*, the selected item will be shown in the viewer.

3.5 Keyword Tagging

In the applications *syngo.share view*, *syngo.share print*, *syngo.share import*, and *syngo.share webview* you can annotate documents. Special authorization from an administrator is needed for this functionality. Tags add another layer of description to documents such that they can be easily found in the future. A configuration determines whether the tags can be chosen from a catalog, dynamically added to a catalog, or entered as free text (e.g. comments). The following DICOM levels support keyword tags: studies, series, and images. Generic containers also support tagging. While searching using keywords, the results are shown independently of the level in which the keyword tag was added.

For further details for tagging read the online help.

3.6 Export/local Saving

In the applications *syngo.share view*, *syngo.share webview*, and *webbox*, there are export functions that allow the data to be saved locally. To export data, the following levels of permission are possible:

- Full permissions export

3.6.1 Exporting with *syngo.share* view

Exporting of studies and series in the document tree:

- Click “Export all series” button to export every file.
- Select the desired items and export them by using the button “Export selected series”.

Exporting of individual images in the print area:

1. Drag the requested element out of the document tree and into the print area.
2. Click the “Export...” button.

3.6.2 Exporting with *syngo.share* webview

Exporting a selection in *syngo.share* webview:

1. Click on either “Export...” or “Export as...” in the context menu.
2. The export wizard will open where you can configure your export.

3.7 Control State of Actions

All desktop applications have a status area, which contain

- If the requested action was successful the state will show a green check mark.
- An exclamation point indicates that the action could not be completed successfully. Click on the status for detailed information.

4 Application *syngo.share* view

The user interface of the application *syngo.share* view is divided into the following sections:

- **Sidebar** on the left
- **Main Area**
- **Document Tree** on the right which can be shown or hidden.

4.1 Display Mode

You can define the following display modes in the lower left sidebar:

- ☰ Full screen display
Display the images in full screen.
- ☰ Multiscreen
To operate with several screens.
- ☰ White screen
To watch a conventional X-ray image all screens except the main screen are getting white.

4.2 Load Area

In this area, you can search for and open patient files either in the archive, the file system, or patient CD. When you select a record in the list of results, all available images for that record will be displayed in the preview window below the list of results. Upon mouse-over, a small magnifying glass icon in the upper lefthand corner of each image preview that will allow you to quickly view images without loading them into the document tree on the righthand side. You can reorder the columns in your search results list by dragging them. You can right-click an item in the results list to get access to a context menu with the following functions:

- ☰ Grouping
To get a better overview, results can be grouped into a maximum of five columns.
- ☰ Search for prior examinations (Patient)
All pre-examinations for each patient are displayed.
- ☰ Search for examinations (Patient)
All examinations for each patient are displayed.
- ☰ Search for prior examinations (Visit)
Examinations are ordered by the number of visits.
- ☰ Acquire using patient ID
Open *syngo.share* import with selected patient ID.
- ☰ Acquire using visit ID
Open *syngo.share* import with selected visit ID.
- ☰ Open
Opens item in the view area and loads documents into the document tree (offers the same functionality as double-clicking the item).
- ☰ Change document properties
Change the description and study data and time of a document.
- ☰ Tag...
Add a keyword tag to a document.
- ☰ Move
You can relocate documents into another organizational unit.
- ☰ Delete
You can delete documents. Depending on the configuration in *syngo.share* view, your selection will either be deleted from the search results only, or from the archive itself.
- ☰ Copy
You can copy documents to another storage area.

Share

You can share documents with colleagues and local physicians, who can view the data with the application *webbox*.

Register Reference Pointer

HL7 ORU messages can be sent to third-party systems in order to inform them about what content is stored in *syngo.share* view.

Change Patient/Visit

Change the patient or visit data for a patient.

Reassign...

This appears instead of Change Patient/Visit if the patient is an entity (usually only used for research purposes).

Send

You can send documents (studies, series, images) to an AET (any DICOM-compatible device or application)

Note

User roles and their associated permissions assigned by the system administrator determine which functions are accessible.

Documents can be loaded into the viewer and document tree by doing either of the following:

- Double-clicking a record in the search results list.
- Selecting a record from the list and clicking the button "Open".

4.3 View Area

In view area, you can view images in a variety of ways depending on the kind of file opened. The following view options are available:

4.3.1 Single-frame view

In the single-frame view, you can browse the images of a series manually by using the toolbar. Every supported file format can be displayed in this mode.

4.3.2 Movie mode

The movie mode provides playback functionality for image series such as play, pause, stop, etc.

4.3.3 Metadata

The metadata view displays meta information of the selected item (e.g. a DICOM image, structured report, etc.). You can have the content of a DICOM/Exif dump displayed as a written representation. You can search for keyword tags, and, for DICOM studies with more than one image, you can select the desired image from the index selection list.

4.3.4 Create and Manage collection folders

You can put images into collections for quick access or in preparation for review. These folders can be created for private use, or you can share them in defined groups. Collection folders only contain references to images. Therefore, when an image is deleted from a collection folder, it will not be deleted from the actual archive. You can create a collection with the following steps:

1. To create a new collection, click "+" in the document tree.
2. Drag the desired images out of the image preview and into the newly created collection folder.
3. Save your data by using the "Save" button, and provide a name, description, and set the visibility of the folder. It is also possible to define the folder as "private".
4. You can search for the folder and make changes by clicking the "collections" button.

4.4 Print Area

In the print area, layout settings can be configured for the printing and exporting of data:

- Click on the print area in the lefthand menu pane.
- Drag the requested items into the print view.
- Click on the button "Print...".

Alternatively you can choose the camera at the context menu at section "view" to display the requested file at the printer view.

It is your decision if you print by a system printer on paper or if you use a special printer for X-ray films. Detailed information you will find at the online help.

EC Declaration of Conformity

Product Details

| | |
|---|------------------|
| Product Name | syngo.share view |
| Type | VA25B |
| Model Number | 11250004 |
| Classification according to Directive 93/42/EEC, Annex IX | Class I |

Manufacturer

| | |
|---------|--|
| Address | ITH icoserve technology for healthcare GmbH Innrain 98 6020 Innsbruck Austria Phone: +43 512 89059 |
|---------|--|

Assessment Details

| | |
|--|-----------|
| Procedure according to Directive 93/42/EEC | Annex VII |
|--|-----------|

Applied Standards

EN ISO 13485:2012 / AC:2012
EN ISO 9001:2008/AC:2009
EN ISO 14971:2012
EN 62304:2006/AC:2008
EN 62366:2008

We declare under sole responsibility that the products described above are in compliance with Directive 93/42/EEC, if their installation and use occurs corresponding to the user and product documentation. This declaration supersedes any declaration previously issued for this product.

Innsbruck, 2018-02-26

Dr. Bernhard Hirsch
Managing Director

Mag. (FH) Herbert Plankl
Head of Quality & Documentation

EC Declaration of Conformity

Manufacturer's Declaration of Conformity
Australian Therapeutic Goods (Medical Devices) Regulations 2002
Full Quality Assurance Procedure

This is a declaration made in accordance with the requirements of Clause 1.8 of Schedule 3 of the **Australian Therapeutic Goods (Medical Devices) Regulations 2002** relating to the stated device.

Reference

syngo.share view

Legal Manufacturer Business Address

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

Facility

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

Medical Device(s)

syngo.share view VA25B
11250004

Authorised Signatory

2018-02-26

2018-02-26

Dr. Bernhard Hirsch
Managing Director

Mag. (FH) Herbert Plankl
Head of Quality & Documentation

Declaration of Conformity AU

- ☐ Classification (per TGA)
 - Class I
 - (Schedule 2, Part 4, Rule 4.1 of the Therapeutic Goods (Medical Devices) Regulations 2002)
- ☐ GMDN Code and Term
 - 41670: Radiology picture archiving and communication system application software
- ☐ Scope of Application
 - This declaration of conformity covers *syngo.share* view:
 - Model Number: 11250004
 - Product Description: *syngo.share* view VA25B

Each kind of medical device to which the Full Quality Assurance Procedures have been applied complies with the applicable provisions of the essential principles, the classification rules, at each stage, from the design of the device until its final inspection before being supplied.

Full Quality Management System Certificate

EN ISO 13485:2012
Certificate No.: 17062224

Standards Applied

(European Harmonised Standard)

- EN ISO 9001:2008, AC:2009
- EN ISO 13485:2012, AC:2012
- EN ISO 14971:2012
- EN 62304:2006, AC:2008
- EN 62366:2008

ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2018

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

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