syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share import

User Manual

Language Directory

Dansk	3
Deutsch	12
Ελληνικά	21
English	30
Español	39
Français	48
Italiano	57
Magyar	66
Nederlands	75
Norsk	84
Polski	93
Português (Brasil)	102
Português (Portugal)	111
Suomi	120
Türkçe	129
简体中文	138

syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share import

Brugermanual

Indhold

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ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2021

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

Made in Austria

syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share import

User Manual

Contents


1	Introduction	32
1.1	Intended Use	32
1.2	Safety Notices	33
	Data Import	33
	Medical Education	33
	Accuracy of Calculations	33
	Units	33
	Access Restriction	33
	Effects of External Programs	34
2	Installation	35
2.1	Hardware Requirements	35
2.2	Software Requirements	35
2.3	Additional System Requirements	35
3	Usage Instructions	36
3.1	Log in with Your User Credentials	36
3.2	Set a Patient Context	36
3.3	Select the Desired Data	36
3.4	View and Edit Data	36
3.5	Select an Archiving Destination	37
3.6	Add Keywords	37
3.7	Complete Archiving	37

1 Introduction

1.1 Intended Use

syngo.share import is software for the electronic conversion and transfer of data and provides the ability to import clinical data from various sources (such as from local filesystems, DICOMDIR media, digital video cameras, digital photo cameras, documents, and DICOM scanners) into connected archives (such as Vendor Neutral Archives, Picture Archiving and Communication Systems, DICOM archives). This product is not intended for diagnosis, prevention, monitoring, prediction, prognosis, treatment, cure or alleviation of diseases or other medical conditions.

1.2 Safety Notices


Data Import

syngo.share import is only intended for the importing of data from various sources.

Medical Education

The use of *syngo.share* import should be granted to qualified staff members only.

Accuracy of Calculations


Images may contain calibrations. Please note that the accuracy of these calibrations correlate to the precision of the program or modality with which they were created. Calibrations must be checked according to the manufacturer's instructions. Calculations, calibrations and results derived therefrom must be checked for plausibility by clinical professionals. Inaccuracies may arise during image capture by clinical professionals or due to the image capturing technology of a given modality. For instance, X-ray images may contain inaccuracies due to the fact that the x-rayed body part is projected onto a two-dimensional imaging surface. Inaccuracies during visualization may be due to hardware settings.

This product rounds calculations as follows: Length and area: if pixel spacing ≥ 0.1 mm: 1 decimal place, if pixel spacing < 0.1 mm: 2 decimal places; density (HU): no decimal place; density average (HU): 1 decimal place; density area (HU): no decimal place. This product uses an interpolation algorithm for image display. This product may not be able to calculate the area of polygons whose edges overlap; this will be indicated by the value "n/a".


Units

The application *syngo.share* import conforms to Directive 80/181/EEC, which contains specifications related to units of calculation. Calibrated lengths are calculated using millimeters. Surface areas are calculated in square millimeters (mm²). Plane angles are calculated in degrees. Densities are calculated using the industry-standard unit, Hounsfield units (HU), which are not regulated by Directive 80/181/EEC.

Access Restriction

Be sure to log out of *syngo.share import* when finished working in order to avoid unauthorized access to the system.

Effects of External Programs

The software manufacturer does not guarantee that third-party applications will not negatively affect the system. The standards of any external programs accessed by *syngo.share import* are set by their respective manufacturers.

2 Installation

Usually, *syngo.share import* is installed via the software deployment mechanism of your hospital. If you install the application yourself, execute the installation file and follow the instructions of the installation guide. For the initial installation, administrator rights are required.

If you need information regarding host, port, and tenant, please contact your system administrator. In case *syngo.share import* is installed without *syngo.share core*, these connection parameters are not required.

2.1 Hardware Requirements

Table 1: Hardware Requirements *syngo.share import*

Part	Requirements
Processor	≥ 2 Cores
Memory	≥ 4 GB
Storage	≥ 80 GB
Network Interface	100BASE-T Ethernet NIC
Recommended Screen Resolution	≥ 1280x1024

2.2 Software Requirements

- Windows 10 Enterprise / Professional Editionx64
- Windows 10 Enterprise / Professional Editionx86

2.3 Additional System Requirements

A paper scanning device must support a TWAIN interface that matches the architecture (32 or 64 bit) of the *syngo.share import* installation. Video sequence and still-image video sources with DirectShow or Media Foundation are supported. Import from DICOM archives with DICOM C-FIND and C-MOVE is supported. DICOM files can be imported from a DICOM-compliant DICOMDIR.

3 Usage Instructions

In the application *syngo.share* import, you can archive any kind of file (DICOM objects, pictures, videos, documents, etc).

Use **F1** to open the Help for more information and a complete list of shortcuts.

To start the archiving process, execute the following steps.

- [Log in with Your User Credentials](#)
- [Set a Patient Context](#)
- [Select the Desired Data](#)
- [View and Edit Data](#)
- [Select an Archiving Destination](#)
- [Add Keywords](#)
- [Complete Archiving](#)

3.1 Log in with Your User Credentials

If necessary, log in with your username and password to use *syngo.share* import.

3.2 Set a Patient Context

Use **F3** to open the patient context tab. You can search for an existing patient or visit, select a patient from a worklist or create a new patient context.

3.3 Select the Desired Data

- In the **Filesystem** tab, you can import data from disk drives and USB flash drives:
 - Double-click on the desired files.
 - Drag them into the loaded documents list.
- You can drag and drop files directly from Windows Explorer.
- To import patient CDs, use the **DICOM Overview** function.
- If a scanner is connected, you can digitalize analog data such as clinical reports or printed X-ray images. They can be edited after scanning and will be added to the archive.
- If a video source is connected, you can record video streams in the **Video** tab. Additionally, you can take a screenshot during a playing videostream by pressing the **F5** key.
- You can load data from any available DICOM node.
- You can also load data from *syngo.share* core.

3.4 View and Edit Data

Use **F10** to switch to view mode. There you can also edit the loaded data, e.g. trim videos and set bookmarks for viewing, or annotate images.

3.5 Select an Archiving Destination

Choose an archiving destination. This can either be a destination in *syngo.share* core or any available DICOM node.

3.6 Add Keywords

Depending on the selected destination, you can assign keywords from certain predefined catalogs or use free text tags.

3.7 Complete Archiving

Use or click on to archive the data.

ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2021

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

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