syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share view

User Manual

Language Directory

Dansk	3
Deutsch	12
Ελληνικά	21
English	31
Español	40
Français	49
Italiano	58
Magyar	67
Nederlands	76
Norsk	85
Polski	94
Português (Brasil)	103
Português (Portugal)	112
Suomi	121
Türkçe	130
简体中文	139

syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share view

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Indhold

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ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2021

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

Made in Austria

syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share view

User Manual English

Contents

1	Introduction	33
1.1	Intended Use	33
1.2	Safety Notices	34
	Image and Result Distribution	34
	Medical Education	34
	Accuracy of Calculations	34
	Units	34
	Access Restriction	34
	Effects of External Programs	35
	Printing of Images	35
	Network Connection	35
	Workstation Accessibility	35
2	Installation	36
2.1	Hardware Requirements	36
2.2	Software Requirements	36
2.3	Additional System Requirements	36
3	Usage Instructions	37
3.1	Log in with Your User Credentials	37
3.2	Search for Patients and Documents	37
3.3	Customize the Results List and Preview the Results	37
3.4	Manage Documents and Metadata	37
3.5	View Documents	38
3.6	Create Collections	38

1 Introduction

1.1 Intended Use

syngo.share view is software to be used in clinical image and result distribution and provides access to clinical data from connected systems (such as Vendor Neutral Archives, Picture Archiving and Communication Systems, DICOM archives, XDS-I Affinity Domains, Electronic Medical Record Systems, Clinical Information Systems) and from local filesystems for the electronic display of data. This product is not intended for diagnosis, prevention, monitoring, prediction, prognosis, treatment, cure or alleviation of diseases or other medical conditions.

1.2 Safety Notices

Image and Result Distribution

syngo.share view is only intended for clinical image and result distribution.
Do not use *syngo.share* view for direct diagnoses.

Medical Education


The use of *syngo.share* view should be granted to qualified staff members only.

Accuracy of Calculations


Images may contain calibrations. Please note that the accuracy of these calibrations correlate to the precision of the program or modality with which they were created. Calibrations must be checked according to the manufacturer's instructions. Calculations, calibrations and results derived therefrom must be checked for plausibility by clinical professionals. Inaccuracies may arise during image capture by clinical professionals or due to the image capturing technology of a given modality. For instance, X-ray images may contain inaccuracies due to the fact that the x-rayed body part is projected onto a two-dimensional imaging surface. Inaccuracies during visualization may be due to hardware settings.

This product rounds calculations as follows: Length and area: if pixel spacing ≥ 0.1 mm: 1 decimal place, if pixel spacing < 0.1 mm: 2 decimal places; density (HU): no decimal place; density average (HU): 1 decimal place; density area (HU): no decimal place. This product uses an interpolation algorithm for image display. This product may not be able to calculate the area of polygons whose edges overlap; this will be indicated by the value "n/a".


Units

The application *syngo.share* view conforms to Directive 80/181/EEC, which contains specifications related to units of calculation. Calibrated lengths are calculated using millimeters. Surface areas are calculated in square millimeters (mm²). Plane angles are calculated in degrees. Densities are calculated using the industry-standard unit, Hounsfield units (HU), which are not regulated by Directive 80/181/EEC.

Access Restriction

Be sure to log out of *syngo.share* view when finished working in order to avoid unauthorized access to the system.

Effects of External Programs

The software manufacturer does not guarantee that third-party applications will not negatively affect the system. The standards of any external programs accessed by *syngo.share* view are set by their respective manufacturers.

Printing of Images

When printing images using the image caption option "As camera", it is possible that some labels will be hidden from view due to the physically limited space on the page. As a result, patient information may not be visible on the printout. Possible transformations such as mirroring, rotation, and applied presentation states are not marked explicitly in the printout. If a series has active presentation states in both the view and print layouts, deactivation of the presentation states in the view layout will also deactivate it in the print layout of the series.

Network Connection

syngo.share view requires an active network connection. For IT security reasons, a (temporary) failure of the network connection may result in you being logged out and having to log in again.

Workstation Accessibility

When using *syngo.share* view, make sure that you can operate the application at all times.

2 Installation

Usually, *syngo.share* view is installed via the software deployment mechanism of your hospital. If you install the application yourself, execute the installation file and follow the instructions of the installation guide. For the initial installation, administrator rights are required.

If you need information regarding host, port, and tenant, please contact your system administrator. In case *syngo.share* view is installed without *syngo.share* core, these connection parameters are not required.

2.1 Hardware Requirements

Table 1: Hardware Requirements *syngo.share* view

Part	Requirements
Processor	≥ 2 Cores
Memory	≥ 4 GB
Storage	≥ 80 GB
Network Interface	100BASE-T Ethernet NIC
Recommended Screen Resolution	≥ 1280x1024

2.2 Software Requirements

- Windows 10 Enterprise / Professional Editionx64
- Windows 10 Enterprise / Professional Editionx86

2.3 Additional System Requirements

Note Regarding Videos

To allow for the highest possible loading speed when viewing videos, the following is recommended:

- The HTTP servers should be activated on FilesystemServers where video data is stored.
- Videos should be imported to an organizational unit for which an archiving rule without data compression is defined.

3 Usage Instructions

syngo.share view is an application that allows you to view and distribute medical documents (e.g. images, reports, videos and photographs). It offers a variety of built-in viewing and annotation tools. You can share, send and print documents. In addition, you can edit and manage patient and document metadata.

Use **F1** to open the Help for more information and a complete list of shortcuts.

3.1 Log in with Your User Credentials

Provided you have not been logged in by a third-party application, use your user credentials to log in to *syngo.share* view.

3.2 Search for Patients and Documents

To begin a search, enter the desired search criteria in the appropriate fields and press the **Search** button. The **Search** button will remain inactive if the entered search criteria are insufficient.

You can customize the visibility of search fields according to your needs and save search queries with predefined criteria as filters. Moreover, you can favorite a filter, which will execute it automatically upon login.

You can search for patients and documents in *syngo.share* core as well as other available DICOM nodes. You can also start a cross-enterprise search (XDS search) for external patients and load external documents, provided that an eHealth Solutions Consumer Adaptor has been configured for this purpose.

3.3 Customize the Results List and Preview the Results

You can arrange the various columns as desired and change the sorting order of the results. You can also group the results according to patients and documents.

Your adjustments are saved in your user profile and will be applied to all future search queries.

Selecting an entry from the results list provides you with detailed series information. Various display options are available.

3.4 Manage Documents and Metadata

Depending on your user permissions, various functions are available for managing documents, patients and their metadata. These functions can usually also be accessed on the series and instance level:

- Open documents in an external application
- Add, edit or delete metadata and keywords
- Export documents (anonymized export is available)
- Send documents to remote peers (e.g. DICOM node)
- Share documents to other applications (e.g. Patient Portal, Inbox, EMR systems)
- Copy documents, move them to another destination or assign them to another patient/visit

- Mark documents as “deleted”

3.5 View Documents

Double-click on a document to open and display it in a predefined layout. You can customize the layout according to your needs.

Different viewing and annotation tools are available depending on the file type:

- Panning, windowing and scrolling functions
- Various transformation and zoom tools
- Various length, area and angle annotation tools
- Drawing and text annotation tools
- Image presets and display Presentation States and DICOM Overlays
- Add, edit or delete metadata and keywords
- Playback tools in *Cine Mode*
- Synchronization of series
- Display of Reference Pointers
- Show document header information

3.6 Create Collections

Use collections to compile documents for personal use or in preparation for review. Collections only contain references to documents, meaning that when you delete a collection, the actual documents will not be deleted from the archive.

You can create private collections or share them with predefined groups of users.

ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2021

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

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