Release VA32C_BF04 / 2024-07-09 / Revision 9813

User Manual

Language Directory

Dansk	3
Deutsch	11
Ελληνικά	19
English	27
Español	35
Français	43
Italiano	51
Nederlands	59
Norsk	67
Polski	75
Português (Portugal)	83
Português (Brasil)	91
Suomi	99
Türkçe	107
简体中文	115

Release VA32C_BF04 / 2024-07-09 / Revision 9813

Brugermanual

Indhold

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User Manual

Contents

1	Intro	oduction	29
	1.1	Intended Use	29
	1.2	Safety Notices	30
		Data Import	30
		Medical Education	30
		Accuracy of Calculations	30
		Units	30
		Access Restriction	30
		Effects of External Programs	30
2	Insta	allation	32
	2.1	Hardware Requirements	32
	2.2	Software Requirements	32
	2.3	Additional System Requirements	32
3	Usag	ge Instructions	33
	3.1	Log in with Your User Credentials	33
	3.2	Set a Patient Context	33
	3.3	Select the Desired Data	33
	3.4	View and Edit Data	33
	3.5	Select an Archiving Destination	34
	3.6	Add Keywords	34
	3.7	Complete Archiving	34

1 Introduction

1.1 Intended Use

syngo.share import is software for the electronic conversion and transfer of data and provides the ability to import clinical data from various sources (such as from local filesystems, DICOMDIR media, digital video cameras, digital photo cameras, documents, and DICOM scanners) into connected archives (such as Vendor Neutral Archives, Picture Archiving and Communication Systems, DICOM archives). This product is not intended for diagnosis, prevention, monitoring, prediction, prognosis, treatment, cure or alleviation of diseases or other medical conditions.

1.2 Safety Notices

Data Import

syngo.share import is only intended for the importing of data from various sources.

Medical Education

The use of syngo.share import should be granted to qualified staff members only.

Accuracy of Calculations

Images may contain calibrations. Please note that the accuracy of these calibrations correlate to the precision of the program or modality with which they were created. Calibrations must be checked according to the manufacturer's instructions. Calculations, calibrations and results derived therefrom must be checked for plausibility by clinical professionals. Inaccuracies may arise during image capture by clinical professionals or due to the image capturing technology of a given modality. For instance, X-ray images may contain inaccuracies due to the fact that the x-rayed body part is projected onto a two-dimensional imaging surface. Inaccuracies during visualization may be due to hardware settings.

This product rounds calculations as follows: Length and area: if pixel spacing \geq 0.1mm: 1 decimal place, if pixel spacing < 0.1mm: 2 decimal places; density (HU): no decimal place; density average (HU): 1 decimal place; density area (HU): no decimal place. This product uses an interpolation algorithm for image display. This product may not be able to calculate the area of polygons whose edges overlap; this will be indicated by the value "n/a".

Units

The application *syngo*.share import conforms to Directive 80/181/EEC, which contains specifications related to units of calculation.

- Calibrated lengths are calculated in millimeters (mm). Uncalibrated lengths are calculated in pixels (pixel).
- Calibrated surface areas are calculated in square millimeters (mm²). Uncalibrated surface areas are calculated in square pixels (pixel²).
- Plane angles are calculated in degrees (°).
- Densities are calculated using the industry-standard Hounsfield unit (HU), which is not regulated by Directive 80/181/EEC.

Access Restriction

Be sure to log out of *syngo*.share import when finished working in order to avoid unauthorized access to the system.

Effects of External Programs

The software manufacturer does not guarantee that third-party applications will not negatively affect the system. The standards of any external programs accessed by *syngo*.share import are set by their respective manufacturers.

2 Installation

Usually, *syngo*.share import is installed via the software deployment mechanism of your hospital. If you install the application yourself, execute the installation file and follow the instructions of the installation guide. For the initial installation, administrator rights are required.

If you need information regarding host, port, and tenant, please contact your system administrator. In case *syngo*.share import is installed without *syngo*.share core, these connection parameters are not required.

2.1 Hardware Requirements

Part	Requirements
Processor	≥ 2 Cores
Memory	≥ 4 GB
Storage	≥ 80 GB
Network Interface	100BASE-T Ethernet NIC
Recommended Screen Resolution	≥ 1280x1024

Table 1: Hardware Requirements syngo.share import

2.2 Software Requirements

One of the following operating systems is required:

- Windows 10 Enterprise / Professional Edition x64 (≥ 1809)
- Windows 11 Enterprise / Professional Edition x64

2.3 Additional System Requirements

A paper scanning device must support a TWAIN interface that matches the 64 bit architecture of the *syngo*.share import installation. Video sequence and still-image video sources with DirectShow or Media Foundation are supported. Import from DICOM archives with DICOM C-FIND and C-MOVE is supported. DICOM files can be imported from a DICOM-compliant DICOMDIR.

3 Usage Instructions

In the application *syngo*.share import, you can archive any kind of file (DICOM objects, pictures, videos, documents, etc).

Use F1 to open the Help for more information and a complete list of shortcuts.

To start the archiving process, execute the following steps.

- Log in with Your User Credentials
- Set a Patient Context
- Select the Desired Data
- View and Edit Data
- Select an Archiving Destination
- Add Keywords
- Complete Archiving

3.1 Log in with Your User Credentials

If necessary, log in with your username and password to use syngo.share import.

3.2 Set a Patient Context

Use **F3** to open the patient context tab. You can search for an existing patient or visit, select a patient from a worklist or create a new patient context.

3.3 Select the Desired Data

- In the Filesystem tab, you can import data from disk drives and USB flash drives:
 - Double-click on the desired files.
 - Drag them into the loaded documents list.
- You can drag and drop files directly from Windows Explorer.
- To import patient CDs, use the **DICOM Overview** function.
- If a scanner is connected, you can digitalize analog data such as clinical reports or printed X-ray images. They can be edited after scanning and will be added to the archive.
- If a video source is connected, you can record video streams in the **Video** tab. Additionally, you can take a screenshot during a playing videostream by pressing the **F5** key.
- You can load data from any available DICOM node.
- You can also load data from syngo.share core.

3.4 View and Edit Data

Use **F10** to switch to view mode. There you can also edit the loaded data, e.g. trim videos and set bookmarks for viewing, or annotate images.

3.5 Select an Archiving Destination

Choose an archiving destination. This can either be a destination in *syngo*.share core or any available DICOM node.

3.6 Add Keywords

Depending on the selected destination, you can assign keywords from certain predefined catalogs or use free text tags.

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İthalatçı ve Türkiye Yetkili Temsilcisi Siemens Healthcare Sağlık A.Ş. Yakacık Caddesi No 111 34870 Kartal İstanbul Türkiye © ITH icoserve technology for healthcare GmbH / Siemens Healthcare GmbH 2024. All rights reserved.
Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen Germany Phone: +49 9131 84-0 siemens-healthineers.com
ITH icoserve technology for healthcare GmbH Innrain 98 6020 Innsbruck Austria Phone: +43 512 89059