**Release VA32C\_BF04 / 2024-07-09 / Revision 9813**

# **syngo.share view**

**User Manual**

# Language Directory

|                      |     |
|----------------------|-----|
| Dansk                | 3   |
| Deutsch              | 11  |
| Ελληνικά             | 19  |
| English              | 27  |
| Español              | 35  |
| Français             | 43  |
| Italiano             | 51  |
| Nederlands           | 59  |
| Norsk                | 67  |
| Polski               | 75  |
| Português (Portugal) | 83  |
| Português (Brasil)   | 91  |
| Suomi                | 99  |
| Türkçe               | 107 |
| 简体中文                 | 115 |

**Release VA32C\_BF04 / 2024-07-09 / Revision 9813**

# **syngo.share view**

**Brugermanual**

# Indhold

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**User Manual**

# Contents

|          |                                                              |           |
|----------|--------------------------------------------------------------|-----------|
| <b>1</b> | <b>Introduction</b>                                          | <b>29</b> |
| 1.1      | Intended Use . . . . .                                       | 29        |
| 1.2      | Safety Notices . . . . .                                     | 30        |
|          | Image and Result Distribution . . . . .                      | 30        |
|          | Medical Education . . . . .                                  | 30        |
|          | Accuracy of Calculations . . . . .                           | 30        |
|          | Units . . . . .                                              | 30        |
|          | Access Restriction . . . . .                                 | 30        |
|          | Effects of External Programs . . . . .                       | 30        |
|          | Printing of Images . . . . .                                 | 31        |
|          | Network Connection . . . . .                                 | 31        |
|          | Workstation Accessibility . . . . .                          | 31        |
| <b>2</b> | <b>Installation</b>                                          | <b>32</b> |
| 2.1      | Hardware Requirements . . . . .                              | 32        |
| 2.2      | Software Requirements . . . . .                              | 32        |
| 2.3      | Additional System Requirements . . . . .                     | 32        |
| <b>3</b> | <b>Usage Instructions</b>                                    | <b>33</b> |
| 3.1      | Log in with Your User Credentials . . . . .                  | 33        |
| 3.2      | Search for Patients and Documents . . . . .                  | 33        |
| 3.3      | Customize the Results List and Preview the Results . . . . . | 33        |
| 3.4      | Manage Documents and Metadata . . . . .                      | 33        |
| 3.5      | View Documents . . . . .                                     | 34        |
| 3.6      | Create Collections . . . . .                                 | 34        |

# 1 Introduction

## 1.1 Intended Use

*syngo.share* view is software to be used in clinical image and result distribution and provides access to clinical data from connected systems (such as Vendor Neutral Archives, Picture Archiving and Communication Systems, DICOM archives, XDS-I Affinity Domains, Electronic Medical Record Systems, Clinical Information Systems) and from local filesystems for the electronic display of data. This product is not intended for diagnosis, prevention, monitoring, prediction, prognosis, treatment, cure or alleviation of diseases or other medical conditions.

## 1.2 Safety Notices

### Image and Result Distribution

*syngo.share* view is only intended for clinical image and result distribution.

**Do not use *syngo.share* view for direct diagnoses.**

### Medical Education

The use of *syngo.share* view should be granted to qualified staff members only.

### Accuracy of Calculations

Images may contain calibrations. Please note that the accuracy of these calibrations correlate to the precision of the program or modality with which they were created. Calibrations must be checked according to the manufacturer's instructions. Calculations, calibrations and results derived therefrom must be checked for plausibility by clinical professionals. Inaccuracies may arise during image capture by clinical professionals or due to the image capturing technology of a given modality. For instance, X-ray images may contain inaccuracies due to the fact that the x-rayed body part is projected onto a two-dimensional imaging surface. Inaccuracies during visualization may be due to hardware settings.

This product rounds calculations as follows: Length and area: if pixel spacing  $\geq 0.1\text{mm}$ : 1 decimal place, if pixel spacing  $< 0.1\text{mm}$ : 2 decimal places; density (HU): no decimal place; density average (HU): 1 decimal place; density area (HU): no decimal place. This product uses an interpolation algorithm for image display. This product may not be able to calculate the area of polygons whose edges overlap; this will be indicated by the value "n/a".

### Units

The application *syngo.share* view conforms to Directive 80/181/EEC, which contains specifications related to units of calculation.

- Calibrated lengths are calculated in millimeters (mm). Uncalibrated lengths are calculated in pixels (pixel).
- Calibrated surface areas are calculated in square millimeters ( $\text{mm}^2$ ). Uncalibrated surface areas are calculated in square pixels ( $\text{pixel}^2$ ).
- Plane angles are calculated in degrees ( $^\circ$ ).
- Densities are calculated using the industry-standard Hounsfield unit (HU), which is not regulated by Directive 80/181/EEC.

### Access Restriction

Be sure to log out of *syngo.share* view when finished working in order to avoid unauthorized access to the system.

### Effects of External Programs

The software manufacturer does not guarantee that third-party applications will not negatively affect the system. The standards of any external programs accessed by *syngo.share* view are set by their respective manufacturers.

### Printing of Images

When printing images using the image caption option “As camera”, it is possible that some labels will be hidden from view due to the physically limited space on the page. As a result, patient information may not be visible on the printout.

Possible transformations such as mirroring, rotation, and applied presentation states are not marked explicitly in the printout. If a series has active presentation states in both the view and print layouts, deactivation of the presentation states in the view layout will also deactivate it in the print layout of the series.

### Network Connection

*syngo.share* view requires an active network connection. For IT security reasons, a (temporary) failure of the network connection may result in you being logged out and having to log in again.

### Workstation Accessibility

When using *syngo.share* view, make sure that you can operate the application at all times.

## 2 Installation

Usually, *syngo.share* view is installed via the software deployment mechanism of your hospital. If you install the application yourself, execute the installation file and follow the instructions of the installation guide. For the initial installation, administrator rights are required.

If you need information regarding host, port, and tenant, please contact your system administrator. In case *syngo.share* view is installed without *syngo.share* core, these connection parameters are not required.

### 2.1 Hardware Requirements

**Table 1:** Hardware Requirements *syngo.share* view

| Part                          | Requirements           |
|-------------------------------|------------------------|
| Processor                     | ≥ 2 Cores              |
| Memory                        | ≥ 4 GB                 |
| Storage                       | ≥ 80 GB                |
| Network Interface             | 100BASE-T Ethernet NIC |
| Recommended Screen Resolution | ≥ 1280x1024            |

### 2.2 Software Requirements

One of the following operating systems is required:

- Windows 10 Enterprise / Professional Edition x64 (≥ 1809)
- Windows 11 Enterprise / Professional Edition x64

### 2.3 Additional System Requirements

#### Note Regarding Videos

To allow for the highest possible loading speed when viewing videos, the following is recommended: The HTTP servers should be activated on FilesystemServers where video data is stored.

## 3 Usage Instructions

*syngo.share* view is an application that allows you to view and distribute medical documents (e.g. images, reports, videos and photographs). It offers a variety of built-in viewing and annotation tools. You can share, send and print documents. In addition, you can edit and manage patient and document metadata.

Use **F1** to open the Help for more information and a complete list of shortcuts.

### 3.1 Log in with Your User Credentials

Provided you have not been logged in by a third-party application, use your user credentials to log in to *syngo.share* view.

### 3.2 Search for Patients and Documents

To begin a search, enter the desired search criteria in the appropriate fields and press the **Search** button. The **Search** button will remain inactive if the entered search criteria are insufficient.

You can customize the visibility of search fields according to your needs and save search queries with predefined criteria as filters. Moreover, you can favorite a filter, which will execute it automatically upon login.

You can search for patients and documents in *syngo.share* core as well as other available DICOM nodes. You can also start a cross-enterprise search (XDS search) for external patients and load external documents, provided that an eHealth Solutions Consumer Adaptor has been configured for this purpose.

### 3.3 Customize the Results List and Preview the Results

You can arrange the various columns as desired and change the sorting order of the results. You can also group the results according to patients and documents.

Your adjustments are saved in your user profile and will be applied to all future search queries.

Selecting an entry from the results list provides you with detailed series information. Various display options are available.

### 3.4 Manage Documents and Metadata

Depending on your user permissions, various functions are available for managing documents, patients and their metadata. These functions can usually also be accessed on the series and instance level:

- Open documents in an external application
- Add, edit or delete metadata and keywords
- Export documents (anonymized export is available)
- Send documents to remote peers (e.g. DICOM node)
- Share documents to other applications (e.g. Patient Portal, Inbox, EMR systems)
- Copy documents, move them to another destination or assign them to another patient/visit
- Mark documents as “deleted”

## 3.5 View Documents

Double-click on a document to open and display it in a predefined layout. You can customize the layout according to your needs.

Different viewing and annotation tools are available depending on the file type:

- Panning, windowing and scrolling functions
- Various transformation and zoom tools
- Various length, area and angle annotation tools
- Drawing and text annotation tools
- Image presets and display Presentation States and DICOM Overlays
- Add, edit or delete metadata and keywords
- Playback tools in Cine Mode
- Synchronization of series
- Display of Reference Pointers
- Show document header information

## 3.6 Create Collections

Use collections to compile documents for personal use or in preparation for review. Collections only contain references to documents, meaning that when you delete a collection, the actual documents will not be deleted from the archive.

You can create private collections or share them with predefined groups of users.

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© ITH icoserve technology for healthcare GmbH / Siemens Healthcare GmbH 2024. All rights reserved.
Siemens Healthcare GmbH
ITH icoserve technology for healthcare GmbH