Release VA32C_BF04 / 2024-07-09 / Revision 9813

syngo.share webview

User Manual

Language Directory

Dansk	3
Deutsch	12
Ελληνικά	21
English	30
Español	39
Français	48
Italiano	57
Nederlands	66
Norsk	75
Polski	84
Português (Portugal)	93
Português (Brasil)	102
Suomi	111
Türkçe	120
简体中文	129

Release VA32C_BF04 / 2024-07-09 / Revision 9813

Brugermanual

Indhold

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User Manual

Contents

1	Intre	aduction	22
•			52
	1.1	Intended Use	32
	1.2	Hardware Requirements	32
	1.3	Software Requirements	32
	1.4	Additional System Requirements	33
	1.5	Safety Notices	34
		Image and Result Distribution	34
		Medical Education	34
		Accuracy of Calculations	34
		Units	34
		Access Restriction	34
		Effects of External Programs	34
		Printing of Images	35
		Network Connection	35
		Workstation Accessibility	35
2	Usag	ge Instructions	36
	2.1	Access the Application and Log in with Your User Credentials	36
	2.2	Search for Patients and Documents	36
	23	Customize the Results List and Preview the Results	36
	2.5	Manage Documents and Metadata	50
	2.4		57
	2.5	View Documents	37
	2.6	Share Documents via the Inbox	37
	2.7	Archive Documents	38

1 Introduction

1.1 Intended Use

syngo.share webview is a web application to be used in clinical image and result distribution and provides access to clinical data from connected systems (such as Vendor Neutral Archives, Picture Archiving and Communication Systems, DICOM archives, XDS/-I Affinity Domains, Electronic Medical Record Systems, Clinical Information Systems) and the local file system for the electronic display of data. This product is not intended for diagnosis, prevention, monitoring, prediction, prognosis, treatment, cure or alleviation of diseases or other medical conditions.

1.2 Hardware Requirements

Table 1: Hardware Requirements syngo.share webview Client

Part	Requirements
Processor	≥ 2 Cores
Memory	\geq 4 GB
Storage	≥ 80 GB
Network Interface	≥ 100 Mbit/s
Network Delay	\leq 100 ms; \leq 60 ms (recommended)
Recommended Screen Resolution	≥ 1280x1024

syngo.share webview is optimized for use on the following mobile devices supporting iPadOS 16:

- iPad (≥ 2020)
- iPad Air (≥ 2019)
- iPad Pro (≥ 2017)

1.3 Software Requirements

The following web browsers are supported:

- Google Chrome \geq 90
- Microsoft Edge \geq 90
- Mozilla Firefox \geq 91
- Mozilla Firefox $ESR \ge 91$
- Apple Safari ≥ 15

Additional requirements for the backend:

• Red Hat Enterprise Linux 8 (\geq 8.6, Podman 4.2.0)

Recommended Setups

The following setups are recommended:

• Google Chrome (latest version) on Windows 10 x64

• Microsoft Edge (latest version) on Windows 10 x64

The operating system Windows 11 x64 is also supported in combination with the latest version of Chrome or Edge.

Note

For the best user experience and optimal performance with the highest security standards, the installation of updates for the web browser and operating system on a regular basis is strongly encouraged.

Mobile Devices

The following mobile web browser is recommended:

• Apple Safari (iPadOS) ≥ 15

Feature restrictions may apply on mobile devices (e.g. export, download and print). Mobile view is required, touch input is recommended. It is further recommended to launch *syngo*.share webview via home screen shortcut (website icon). This will cause the browser to open in full-screen mode, which offers the best user experience of *syngo*.share webview. How to set a shortcut on the home screen is described in the online iPad user guide.

1.4 Additional System Requirements

Note Regarding Videos

To allow for the highest possible loading speed when viewing videos, the following is recommended:

- The HTTP servers should be activated on FilesystemServers where video data is stored.
- The video codec H.264 (Level: 3.0, Profile: Baseline or Constrained Baseline) and audio codec AAC should be used.

Note Regarding Whole Slide Images for Digital Pathology

DICOM Whole Slide Images for Digital Pathology are only supported in conjunction with *syngo*.share core.

1.5 Safety Notices

Image and Result Distribution

syngo.share webview is only intended for clinical image and result distribution. **Do not use** *syngo*.share webview for direct diagnoses.

Medical Education

The use of syngo.share webview should be granted to qualified staff members only.

Accuracy of Calculations

Images may contain calibrations. Please note that the accuracy of these calibrations correlate to the precision of the program or modality with which they were created. Calibrations must be checked according to the manufacturer's instructions. Calculations, calibrations and results derived therefrom must be checked for plausibility by clinical professionals. Inaccuracies may arise during image capture by clinical professionals or due to the image capturing technology of a given modality. For instance, X-ray images may contain inaccuracies due to the fact that the x-rayed body part is projected onto a two-dimensional imaging surface. Inaccuracies during visualization may be due to browser or hardware settings.

This product rounds calculations as follows: Length and area: 1 decimal place; density (HU): no decimal place; area density (HU): 1 decimal place. This product uses an interpolation algorithm for image display.

Units

The application *syngo*.share webview conforms to Directive 80/181/EEC, which contains specifications related to units of calculation.

- Calibrated lengths are calculated in millimeters (mm). Uncalibrated lengths are calculated in pixels (pixel).
- Calibrated surface areas are calculated in square millimeters (mm²). Uncalibrated surface areas are calculated in square pixels (pixel²).
- For DICOM Whole Slide Images (WSI), the unit display changes to micrometers (μm) / square micrometers (μm^2) if the length falls below 1 mm or the surface area falls below 0.1 mm².
- Plane angles are calculated in degrees (°).
- Densities are calculated using the industry-standard Hounsfield unit (HU), which is not regulated by Directive 80/181/EEC.

Access Restriction

Be sure to log out of *syngo*.share webview when finished working in order to avoid unauthorized access to the system.

Effects of External Programs

The software manufacturer does not guarantee that third-party applications will not negatively affect the system. The standards of any external programs accessed by *syngo*.share webview are set by their respective manufacturers.

Printing of Images

To print, *syngo*.share webview transfers the data to the print function of the browser. Thus there may be differences between the view in the Image Display Area in *syngo*.share webview and the printout, e.g. the image could be cut off, annotations could be incomplete or missing entirely.

Possible transformations such as mirroring, rotation, and applied presentation states are not marked explicitly in the printout.

Network Connection

syngo.share webview requires an active network connection. For IT security reasons, a (temporary) failure of the network connection may result in you being logged out and having to log in again.

Workstation Accessibility

When using *syngo*.share webview, make sure that you can operate the application at all times.

2 Usage Instructions

syngo.share webview is a web application that allows you to view and distribute medical documents (e.g. images, reports, videos and photographs). It offers a variety of built-in viewing and annotation tools. You can upload documents directly to a patient context, as well as share and send them. In addition, you can edit and manage patient and document metadata.

2.1 Access the Application and Log in with Your User Credentials

Provided you have not been logged in by a third-party application, use the following URL template to access *syngo*.share webview:

https://<server>/webview/<tenant>/<site>

Enter your user credentials to log in to *syngo*.share webview. If applicable, you also have to select the domain. If you want to preselect the domain, you can provide this information in the URL according to the following template:

https://<server>/webview/<tenant>/<site>/?domain=<domain>

If your institution only has a single tenant and/or site, you do not have to provide <tenant> and/or <site>. If you are using *syngo*.share webview without *syngo*.share core, tenant and site do not apply.

In case trusted-only login is activated for your institution, access to *syngo*.share webview is only possible via a valid SAML token or HMAC signature.

2.2 Search for Patients and Documents

To begin a search, enter the desired search criteria in the appropriate fields and press the **Search** button. The **Search** button will remain inactive if the entered search criteria are insufficient.

You can customize the visibility of search fields according to your needs and save search queries with predefined criteria as filters. Moreover, you can favorite a filter, which will execute it automatically upon login.

You can search for patients and documents in *syngo*.share core as well as other available DICOM nodes. You can also start a cross-enterprise search (XDS search) for and open external documents for a given patient, provided that an eHealth Solutions Consumer Adaptor has been configured for this purpose.

2.3 Customize the Results List and Preview the Results

You can arrange the various columns as desired and change the sorting order of the results. You can also group the results according to patients and documents as well as filter them via text entry.

Your adjustments are saved in your user profile and will be applied to all future search queries.

Selecting an entry from the results list provides you with detailed series information and referenced documents as well as an interactive preview. Various display options are available.
In addition, you can view available prior examinations for found DICOM document entries and preview their contents. Opening a prior examination in the Viewer displays it side by side with the corresponding current examination.

2.4 Manage Documents and Metadata

Depending on your user permissions, various functions are available for managing documents, patients and their metadata. These functions can usually also be accessed on the series and instance level:

- Open documents in an external application
- Add, edit or delete metadata and keywords
- Export documents (anonymized export is available)
- Send documents to remote peers (e.g. DICOM node)
- Share documents to other applications (e.g. Patient Portal, Inbox, EMR systems)
- Copy documents, move them to another destination or assign them to another patient/visit
- Mark documents as "deleted"
- Rearrange studies and series

2.5 View Documents

Double-click on a document to open and display it in a predefined layout. You can customize the layout according to your needs.

Different viewing and annotation tools are available depending on the file type:

- Panning, windowing and scrolling functions
- Various transformation and zoom tools
- Various length, area and angle annotation tools
- Drawing and text annotation tools
- Image presets and display Presentation States and DICOM Overlays
- Add, edit or delete metadata and keywords
- Playback tools in Cine Mode
- Synchronization of series
- Display of Reference Pointers

The patient banner is always visible in the Viewer to prevent patient mix-ups.

The *Patient Jacket* displays all available documents in the archive for the current patient. You can also search for corresponding documents in third-party archives.

2.6 Share Documents via the Inbox

The inbox is used for the electronic transmission of documents and is available to all registered doctors.

You can share documents with one or more selected users. The recipient(s) are automatically informed when documents are shared with them.

The inbox offers the same functions as the results list (see Section 2.3) and the Viewer (see Section 2.5).

2.7 Archive Documents

syngo.share webview allows you to archive documents directly to patients or visits.

Archiving documents is only supported for *syngo*.share core.

To archive documents, follow these steps:

- **1.** Set a patient context.
- 2. Select the required data.
- 3. Select an archiving destination.
- 4. Add keywords.
- 5. Complete archiving.

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İthalatçı ve Türkiye Yetkili Temsilcisi Siemens Healthcare Sağlık A.Ş. Yakacık Caddesi No 111 34870 Kartal İstanbul Türkiye © ITH icoserve technology for healthcare GmbH / Siemens Healthcare GmbH 2024. All rights reserved.
Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen Germany Phone: +49 9131 84-0 siemens-healthineers.com
ITH icoserve technology for healthcare GmbH Innrain 98 6020 Innsbruck Austria Phone: +43 512 89059