

# Fluorescence lifetime imaging and spectroscopic analysis combined in a surgical microscope for optimized fluorescence-guidance in brain tumor resection

David Reichert<sup>a,b</sup>, Mikael T. Erkkilae<sup>a</sup>, Georg Widhalm<sup>c</sup>, Barbara Kiesel<sup>c</sup>, Lisa Wadiura<sup>c</sup>, Petra Mercea<sup>c</sup>, Thomas Roetzer<sup>d</sup>, Adelheid Woehrer<sup>d</sup>, Marco Andreana<sup>a</sup>, Angelika Unterhuber<sup>a</sup>, Marco Wilzbach<sup>e</sup>, Christoph Hauger<sup>e</sup>, Wolfgang Drexler<sup>a</sup>, and Rainer A. Leitgeb<sup>a,b</sup>

<sup>a</sup>Center for Medical Physics and Biomedical Engineering, Medical University of Vienna, Austria

<sup>b</sup>Christian Doppler Laboratory OPTRAMED, Medical University of Vienna, Austria

<sup>c</sup>Department of Neurosurgery, General Hospital and Medical University of Vienna, Austria

<sup>d</sup>Division of Neuropathology and Neurochemistry, Department of Neurology, Medical University of Vienna, Austria

<sup>e</sup>Carl Zeiss Meditec AG, Advanced Development Microsurgery, Germany

## Objective

- About 25,000 malignant brain tumors are prognosed for 2020 in the United States, ~ 50 % of those glioblastoma with a median survival of ~ 8 months<sup>1</sup>
- Extent of resection is a key prognostic factor in the treatment of brain tumors
- Protoporphyrin IX (PpIX) fluorescence guidance intraoperatively visualizes tumorous tissue and has shown to increase the extent of resection<sup>2</sup>




Figure 1: White-light (left) and fluorescence image (right) of a conventional surgical microscope. Courtesy: Georg Widhalm (AKH Wien)

- Weak PpIX concentrations, as found e.g. in low-grade glioma (LGG) or infiltration zones of high-grade glioma (HGG), cannot be visualized due to spectrally overlapping autofluorescence




Figure 2: Emission spectra of PpIX  
Left: strong PpIX fluorescence; Right: weak PpIX fluorescence on an autofluorescence background  
 $\tau$  = fluorescence lifetime

Objective → Sensitive PpIX detection through a surgical microscope by measuring the fluorescence lifetime  $\tau$  for improved brain tumor resection

## Materials & Methods




Figure 3: Laboratory setup, BF: bandpass filter, 590 – 740 nm; PMT: photomultiplier tube

- Raster-scanning frequency domain fluorescence lifetime imaging was integrated in a surgical microscope together with a spectrometer
- FLIM acquisition time about 16 seconds at a spatial resolution of 40  $\mu$ m
- Spectroscopic information can be acquired on selected regions of interest (ROI)
- Brain tumor samples were imaged ex vivo within 1 hour after resection

## Funding

This project has received funding from the Austrian Christian Doppler Research Association as well as from the innovation board of the Carl Zeiss Meditec AG.

## Results

Figure 4: High-grade glioma - strong visible fluorescence during surgery

- Fluorescence lifetime increased up to 16 ns
- Typical peaks of PpIX at 635 and 705 nm in spectra




Figure 5: Low-grade glioma - no visible fluorescence during surgery

- Weak PpIX signal on an autofluorescence background (ROI 1 & 2) increased the fluorescence lifetime up to 5 ns and could be contrasted from an area with pure autofluorescence ( $\tau < 2$  ns, ROI 3)






Figure 6: Violin plots with the lifetime distribution for LGG & HGG (INF: infiltration zone, REA: reactive tissue, TUM: core tumor)

Alanegativ: no visible fluorescence during surgery  
Alapositiv: visible fluorescence during surgery

Median and quartiles are indicated by solid and dashed lines

Figure 7: PpIX vs. autofluorescence signal contribution at 635 nm plotted as a function of the fluorescence lifetime measured at the respective ROI. Higher PpIX signal contributions entailed increased lifetimes when measured with frequency domain FLIM.



## Conclusion

FLIM enables the delineation of weak PpIX signals from an autofluorescence background. Integrated into a surgical microscope, this technology potentially improves the extent of resection and therefore patient outcome.

## References

1. Ostrom, Q. T. et al. CBTRUS Statistical Report: Primary Brain and Other Central Nervous System Tumors Diagnosed in the United States in 2013–2017. Neuro-Oncol. 22, iv1–iv96 (2020).
2. Stepp, H. & Stummer, W. 5-ALA in the management of malignant glioma. Lasers Surg. Med. 50, 399–419 (2018).
3. Russell, J. A. et al. Characterization of Fluorescence Lifetime of Photofrin and Delta-Aminolevulinic Acid Induced Protoporphyrin IX in Living Cells Using Single- and Two-Photon Excitation. IEEE J. Sel. Top. Quantum Electron. 14, 158–166 (2008).
4. Kattelhardt, S. R. et al. Multiphoton excitation fluorescence microscopy of 5-aminolevulinic acid induced fluorescence in experimental gliomas. Lasers Surg. Med. 40, 273–281 (2008).